



GRAFTING & TECHNOLOGY

CLIP FOR CUCURBITACEAE



**TRIA20M
CLIP**



**TRIA21A
CLIP**



**TWEEZER25C1
CLIP D.25mm CIR/AZZ**



**TWEEZER25M1
CLIP D. 25mm**



**TWEEZER25R1
CLIP D.2,5mm RIG/VE**



**TWEEZER27R
CLIP D.27mm RIG/NEU**



**TWEEZER30
CLIP D.30MM**



INSTRUCTIONS FOR USE:

Rootstock: pumpkin (shintoza, macis, vita)

Scion: watermelon (TWEEZER27R)
melon
cucumber

Graft: with pumpkin cutting, leaving only one cotyledon
cucumber scion needs to be mounted clashed
watermelon with TWEEZER27R: rootstock direct seeding leaving two cotyledons with central cut

Material: non toxic rawmaterial, colored with food coloring